

HTM4100 15" TFT Touch Screen display

TECHNICAL INFORMATION

Medicare display screens are simple, aesthetic and feature a user-friendly touch screen interface as standard.

medicare
SYSTEMS



Display

15 inch colour backlit screen
HTM compliant screen colours
Large font clear text display
Wide angle viewing
Option of numeric or alphanumeric display
(e.g. room numbers or names)

Location maps

Optional map display to show the location of each call on the relevant floor plan
Colour coded location markers flash for call types

Sounder

HTM compliant sounds as standard
Client adjustable day/night volumes,
pre-programmable times

Call history

Record of all calls and events, time and date stamped
Response time *(from call to attendance)*
Attendance time *(nurse presence to cancel)*
Staff identity *(when used with ID fobs)*

Zoning

Fully Zoneable
Text addresses for zones
Zone priority option
Call escalation across zones

Power requirements

230V non-switched socket (or non-switched fused spur)

User-editable information includes

Patient/resident notes
Staff reminders for night checks, medication etc.
Resident/room names
Switchable day/night mode
Assistive technology interface
(for effective dementia management)
Optional eLogging available
(subject to clients broadband)

Call accept

Touch screen activated call accept
The accept feature shows that a member of staff has selected the call to attend

Function buttons

Control buttons on touch screen
Volume adjust up/down/mute
Calls Waiting – to scroll through waiting call list*
Map – to display call location
Notes – display resident notes
Menu – to access advance functions

**Auto scroll option available.
The display will cycle through the first four calls waiting*

Dimensions (mm)

W x D x H: 392 x 80 x 325 (440 with aerial)

Compliance

R & TTE, EMC
Low voltage safety directives

